

One-stop PURCHASE



Professional SER VICE



Paints or inks often are in different shear rates status from manufacturing to application. As a typical Non-newtonian fluid, paints or inks always show different viscosity characteristics at different shear rates.

Generally speaking, paints are in a low shear rate state when it is stored, transported, levelling and sagging. But they are in the medium shear rate when being pumped, dipped and mixed at a low speed. And at high-speed dispersion, roll, spray and brush coating, they are in a high shear state, and the shear rate is generally 9000S⁻¹ to 12000S⁻¹. Therefore, understanding the rheological characteristics of paints or inks at such high shear rates requires the use of cone-plate viscometer to measure.

BGD 182 Cone-and-plate Viscometer uses a special angle conical spindle, which generates a very high shear rate on the measured sample under the high-speed drive of a stepper motor. It conforms to ISO 2884-1 《 Paints and varnishes — Determination of viscosity using rotary viscometers — Part 1: Cone-and-plate viscometer operated at a high rate of shear》.

Features







Features

- ♦ 7-inch touch screen, menu operation, rich display content (measurement values, spindle No., speed, shear rate, etc.), simple and convenient to operate.
- Beautiful and durable metal shell, Use a handle to lift up and down, precise and fast positioning, high reliability.
- The electrical part of the whole machine has been carefully designed, with high reliability, precision, stability, and ease of use.
- ARM chip processor, with faster data processing speed.
- ♦ High precision linear calibration before leaving the factory, capable of interpolating multiple points to ensure higher measurement accuracy.
- ◆ Measure the sample viscosity by stepless speed, convenient for operators to choose different shear rates based on the actual application conditions of the sample.
- The instrument automatically alarms when the measurement exceeds the range.



One-stop PURCHASE



Professional SER VICE

- ♦ The instrument automatically alarms when the measurement exceeds the range.
- ♦ High speed data transmission interface ensures a fast and stable communication connection between the instrument and the computer.
- ◆ Data storage and exporting data by an external USB function.
- Built in high-precision PT100 temperature sensor, with high temperature control accuracy and good stability.
- With calibration function to temperature, ensures precise and reliable temperature for controlling.
- ♦ Powerful and simple calibration function: Users only need one bottle of standard oil and can quickly calibrate the instrument through its built-in calibration menu.
- ◆ Can choose the data collection and programmed analysis software to analyze the rheological features of the sample thoroughly.

Main Technical Parameters:

- ★ Adjustable Speed Range: 5rpm ~ 1000rpm (continuously viable, @1rpm)
- ★ Measurement error: < 2% of full scale
- ★ Sample Volume: < 2ml (See the table below for details)
- ★ Two built-in temperature control for option: 5°C-75°C (L type / Low temperature) ;

 50°C –235 $^{\circ}\text{C}$ (H type / High temperature)

- ★ Temperature resolution: 0.1°C
- ★ Temperature control accuracy: ± 0.5 °C (L type); ± 1 °C (H type and the set temperature is less than 150 °C) or ± 2 °C

(H type and the set temperature is higher than 150° C)

- ★ Rotors: 10 types for option (corresponding shear rates and measurement ranges are shown in the table below)
- ★ Power Supply: AC 220V, 50Hz/60Hz, the max. current is about 1.5A
- ★ Size: 275mm × 210mm × 460mm (L×W×H)
- ★ Net weight: 12KG
- ★ Ordering information: BGD 182/1---Low torque Cone-and-plate Viscometer;

BGD 182/2---High torque Cone-and-plate Viscometer

Spindle type	CAP-01	CAP-02	CAP-03	CAP-04	CAP-05	CAP-06	CAP-07	CAP-08	CAP-09	CAP-10
Sample size	67µL	38µL	24µL	134µL	67μL	30µL	1,700µL	400µL	100µL	170µL
Shear rate range (S ⁻¹)	66.5~ 13,300	66.5~ 13,300	66.5~ 13,300	16.5~ 3,300	16.5~ 3,300	16.5~ 3,300	13~ 2,000	13~ 2,000	13~ 2,000	25~ 5,000
BGD 182/1 Measurement range (mPa.s)	20~ 1,600	20~ 3,200	20~ 6,600	20~ 13,000	20~ 26,000	20~ 66,000	20~ 2,600	20~ 10,800	20~ 44,000	20~ 4,400
BGD 182/2 Measurement range (mPa.s)	20~ 37,500	40~ 75,000	80~ 150,000	100~ 300,000	300~ 600,000	800~ 1,500,000	78~ 62,500	313~ 250,000	125~ 1,000,000	100~ 100,000

Note: Calculating the shear rate: $13.33 \times$ the current speed (For CAP-01 to CAP-03);

 $3.33 \times$ the current speed $\,$ (For CAP-04 to CAP-06) ; 2 \times the current speed $\,$ (For CAP-07 to CAP-09) ;

 $5 \times \text{the current speed (For CAP-10)}$

Tel: +86-20-82169666 E-mail: sale: jack@biuged.com Service: service@biuged.com WhatsApp: +86-13316053991 Website: www.biuged.com